

KLT picker 2.0 – user guide

Amitay Eldar*, Boris Landa, Yoel Shkolnisky†

January, 2022

1 Downloading and running the KLT picker- MATLAB package (source code)

Running the KLT picker from its source code requires MATLAB to be installed on your system.

Download and extract the package from

https://github.com/ShkolniskyLab/kltpicker2_matlab

The extracted directory has the following structure:

KLTpicker_start	Main script
userManual	Documentation for the package
LICENSE	Package license
README	A short README file
matlab/	MATLAB source code
get10028	Script to download micrographs for the example of Section 3.
example/	Directory used in Section 3 to demonstrate the KLT picker

After downloading the package, start MATLAB in the directory of the KLTpicker package and run

`KLTpicker_start`

Using the KLT picker once started is described in Section 2.

2 Using the KLT picker

The following questions will appear, one by one

Enter full path of micrographs MRC file:

Type the path to the directory which contains the micrographs MRC files.

*amitayeldar@tauex.tau.ac.il

†yoelsh@tauex.tau.ac.il

Enter full path of output directory:

Type the path where the particles coordinate files will be saved.

Enter the particle size in pixels:

Type the particle diameter in pixels (more precisely, the diameter of the extracted box).

pick all particles?(Y/N)? [Y]

Type Y to pick all particles using the optimal threshold derived on the paper. If you type N, then the following question will appear

How many particles to pick:

Type the number of particles to pick in each micrograph.

Pick noise images?(Y/N)? [N]

Type N if you don't want to pick noise images. If you type Y then the following question will appear

How many noise images to pick:

Type the number of noise images to pick in each micrograph.

Do you want to use ASOCCEM for contamination removal?(Y/N)? [Y]

Type N if you don't want to use ASOCCEM for contamination removal. If you type Y then the following question will appear

Do you want to save ASOCCEM masks?(Y/N)? [N]

Type Y if you want to save for each micrograph a binary mask (.mrc) which indicates the contaminated regions.

Do you want to change ASOCCEM default parameters?(Y/N)? [N]

If you type Y then the following question will appear

Enter ASOCCEM downsampling image size (should be a positive number):

Type the downsampling image size (the default parameter is 600).

Enter ASOCCEM covariance area size (should be a positive odd number)

Type the covariance area size (the default parameter is 5).

Do you want to use the GPU?(Y/N)? [Y]

Type Y to use the GPUs found on your system, and N to not use them (CPU only).

The KLT picker will then start and will display progress notifications. The outputs are the coordinate files (box and star), a text file summarizing the picking process and the contamination masks (if the user chose to save them) at the output directory.

3 Example

In this section we demonstrate the use of the KLT picker on 5 micrographs of the EMPIAR-10028 data set (Plasmodium Falciparum 80S ribosome) [2] from the EMPIAR repository [1].

If MATLAB is installed, start MATLAB in the package directory and run

```
get10028 # downloading 5 micrographs to ./example/micrographs
directory.
```

If MATLAB is not installed, download the micrographs manually from

<https://www.ebi.ac.uk/pdbe/emdb/empiar/entry/10028/>

Start the KLT picker as described in Section 1. For example, if MATLAB is installed, type

```
KLTpicker_start
```

Once starting the KLT picker, enter the following input

```
Enter full path of micrographs MRC file: ./example/micrographs/
```

```
Enter full path of output directory: ./example/results
```

```
Output directory does not exist. Create?(Y/N)? [Y] Y
```

```
Enter the particle size in pixels: 300
```

```
pick all particles?(Y/N)? [Y] Y
```

```
Pick noise images?(Y/N)? [N] N
```

```
Do you want to use ASOCCEM for contamination removal?(Y/N)? [Y] Y
```

```
Do you want to save ASOCCEM masks?(Y/N)? [N] N
```

```
Do you want to change ASOCCEM default parameters?(Y/N)? [N] N
```

```
Do you want to use the GPU?(Y/N)? [Y] Y
```

Once the KLT picker has finished, it will display the message: **Finished the picking successfully.** In order to display the picking results in EMAN [3], open a terminal in the package directory and create a new directory named **eman**, using the command

```
mkdir eman
```

Change directory to **eman** and enter the following commands one by one:

```

e2rawdata.py ../example/micrographs/*.mrc --invert --edgenorm
--xraypixel --ctfest --apix=1.0 --voltage=200.0 --cs=2.0 --ac=10.0
--threads=4 --defocusmin=0.6 --defocusmax=4.0

e2import.py
../example/results/pickedParticlesParticleSize300/box/*.box
--import_boxes --box_type=boxes

e2boxer.py --allmicrographs --boxsize=300 --ptclsize=300 --apix=1.0
--no_ctf --gui --threads=4

```

In order to display the picking results of the **first** micrograph in RELION [4], open the terminal in the package directory and create a new directory named **relion**, using the command

```
mkdir relion
```

Change directory to **relion** and enter the following command

```

relion_display --i ../example/micrographs/001.mrc --coords
../example/results/pickedParticlesParticleSize300/star/001.star
--scale 0.15 --particle_radius 150 --angpix 1 --lowpass 20 --pick

```

In order to display the picking results of the n'th micrograph change 001.mrc and 001.star in the above command to 00n.mrc and 00n.star.

4 Citation

If you use the KLT picker, please cite *KLT picker: Particle picking using data-driven optimal templates*, *Journal of Structural Biology*, Accepted for publication.. A preprint is available at <https://arxiv.org/abs/1912.06500>.

References

- [1] Iudin, A., Korir, P., Salavert-Torres, J., Kleywegt, G., and Patwardhan, A. (2016). *EMPIAR: A public archive for raw electron microscopy image data*. Nature Methods, 13.
- [2] Wong, Wilson and Bai, Xiao-chen and Brown, Alan and Fernandez, Israel S and Hanssen, Eric and Condrón, Melanie and Tan, Yan Hong and Baum, Jake and Scheres, Sjors H W. (2014). *Cryo-EM structure of the Plasmodium falciparum 80S ribosome bound to the anti-protozoan drug emetine*. Wong et al. eLife ,3,e03080.
- [3] G. Tang, L. Peng, P.R. Baldwin, D.S. Mann, W. Jiang, I. Rees & S.J. Ludtke. (2007). *EMAN2: an extensible image processing suite for electron microscopy*. Journal of structural biology, 157, 38-46.
- [4] Scheres, Sjors HW. (2015). *Semi-automated selection of cryo-EM particles in RELION 1.3*. Journal of structural biology, 128, 114-122.